


OMEGA Thomson Scattering System Upgrade



J. KATZ, R. BONI, D. H. FROULA, G. GATES, A. NAUSS, J. SZCZEPANSKI, M. J. SHOUP III, AND A. AGLIATA

University of Rochester, Laboratory for Laser Energetics

Summary


The Upgraded Thomson-Scattering System is a Robust Diagnostic Platform on the OMEGA Laser



- The Thomson-scattering system measures plasma conditions such as ion and electron temperature, electron density, and plasma-flow velocity
- Fully reflective optical transport provides high-quality imaging for wavelengths ranging from 200 nm to 800 nm
- Improved UV signal throughput allows for the detection of Thomson scattering from electron-plasma wave features using a 263-nm probe beam
- Two streak-camera-coupled spectrometers and a gated optical imager are available to support a variety of experimental objectives


E21091

A New TIM-Based Reflective Telescope Provides High-Quality Imaging of Optical Emission from 200 to 800 nm




E20931

The TIM-Based Telescope Uses an Unobstructed Section of a Schwarzschild Objective




E20932

Initial Telescope Performance Tests Demonstrate High Spatial Resolution and Diffraction-Limited Spot Sizes



E20933

A New Instrument Platform Houses Diagnostics for Time-Resolved Spectroscopy and Gated 2-D Imaging




- Ion-acoustic wave system
 - high-resolution spectrometer coupled to ROSS streak camera
- Electron-plasma wave system
 - broadband spectrometer coupled to ROSS streak camera
- Optical imager
 - high-resolution gated charge-coupled-device (CCD) camera

Each of the three detector systems can be independently configured to support a variety of experimental objectives.

E20934


The Ion-Acoustic Wave System Consists of a 1-m Spectrometer Coupled to a ROSS Streak Camera



Spectrometer focal length	1.0 m
Spectral resolution	0.03 to 0.1 nm
Spectral window	4 to 6 nm
Detectable wavelengths	200 to 800 nm
Temporal resolution	200 ps
Sweep speeds	1.5 ns, 5 ns, 15 ns, 25 ns
Maximum field of view at TCC	370 μm
Image magnification	1.4 \times

E20934


The Electron-Plasma Wave System Consists of a 0.3-m Spectrometer Coupled to a ROSS Streak Camera



Spectrometer focal length	0.3 m
Spectral resolution	0.5 to 2.0 nm
Spectral window	75 to 375 nm
Detectable wavelengths	200 to 800 nm
Temporal resolution	100 ps
Sweep speeds	1.5 ns, 5 ns, 15 ns, 25 ns
Maximum field of view at TCC	275 μm
Image magnification	2.1 \times

E20935

The Optical Imager Consists of a Gated Charge-Coupled-Device Camera



Detector	PI-Max-3 Gated CCD
Detectable wavelengths	200 to 800 nm
Minimum gate duration	3 ns
Field of view at TCC	1.5 mm
Spatial resolution at TCC	20 μm
Image magnification	3.2 \times

E20936


Summary

The Upgraded Thomson–Scattering System is a Robust Diagnostic Platform on the OMEGA Laser





- The Thomson–scattering system measures plasma conditions such as ion and electron temperature, electron density, and plasma-flow velocity
- Fully reflective optical transport provides high-quality imaging for wavelengths ranging from 200 nm to 800 nm
- Improved UV signal throughput allows for the detection of Thomson scattering from electron-plasma wave features using a 263-nm probe beam
- Two streak-camera–coupled spectrometers and a gated optical imager are available to support a variety of experimental objectives

A New TIM-Based Reflective Telescope Provides High-Quality Imaging of Optical Emission from 200 to 800 nm



The TIM-Based Telescope Uses an Unobstructed Section of a Schwarzschild Objective




- Diffraction limited
- $f/10$ collection
- Unobstructed
- Spherical optics
- All-reflective transport
- Aluminum coatings
- $>80\%$ reflectivity 200 to 800 nm

Initial Telescope Performance Tests Demonstrate High Spatial Resolution and Diffraction-Limited Spot Sizes




Imaged microscopy mesh grid



Hole width: $28\ \mu\text{m}$, bar width: $23\ \mu\text{m}$
Image magnification: $2.3\times$
Illumination wavelength: $635\ \text{nm}$


Airy disk formed imaging
 $4\text{-}\mu\text{m}$ core single-mode fiber

Image




Wavelength: $635\ \text{nm}$
Working $f/\#$: 23
Calculated airy disk: $35.6\ \mu\text{m}$

Lineout




A New Instrument Platform Houses Diagnostics for Time-Resolved Spectroscopy and Gated 2-D Imaging



- Ion-acoustic wave system
 - high-resolution spectrometer coupled to ROSS streak camera
- Electron-plasma wave system
 - broadband spectrometer coupled to ROSS streak camera
- Optical imager
 - high-resolution gated charge-coupled-device (CCD) camera


Each of the three detector systems can be independently configured to support a variety of experimental objectives.

The Ion-Acoustic Wave System Consists of a 1-m Spectrometer Coupled to a ROSS Streak Camera




Spectrometer focal length	1.0 m
Spectral resolution	0.03 to 0.1 nm
Spectral window	4 to 6 nm
Detectable wavelengths	200 to 800 nm
Temporal resolution	200 ps
Sweep speeds	1.5 ns, 5 ns, 15 ns, 25 ns
Maximum field of view at TCC	370 μm
Image magnification	1.4\times

Electron-Plasma Wave System Consists of a 0.3-m Spectrometer Coupled to a ROSS Streak Camera



Spectrometer focal length	0.3 m
Spectral resolution	0.5 to 2.0 nm
Spectral window	75 to 375 nm
Detectable wavelengths	200 to 800 nm
Temporal resolution	100 ps
Sweep speeds	1.5 ns, 5 ns, 15 ns, 25 ns
Maximum field of view at TCC	275 μm
Image magnification	2.1\times

The Optical Imager Consists of a Gated Charge-Coupled-Device Camera



Detector	Pi-Max-3 Gated CCD
Detectable wavelengths	200 to 800 nm
Minimum gate duration	3 ns
Field of view at TCC	1.5 mm
Spatial resolution at TCC	20 μm
Image magnification	3.2\times