



X-ray diffraction on Z





Exceptional service in the

national

interest

Tommy Ao and Marius Schollmeier

J. Benage, E. Field, M.W. Kimmel, J.L. Porter, P.K. Rambo, J. Schwarz, and C.T. Seagle

Sandia National Laboratories, Albuquerque, NM

National ICF Diagnostics Working Group Meeting Los Alamos National Laboratory, Los Alamos, NM, October 6-8, 2015





Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.

Combining x-ray diffraction with Z's unique high energy density samples will provide benchmark quality data

- Z's high energy density matter samples are large, uniform, long-lived and precisely characterized
- X-ray diffraction will expand diagnostic capabilities on Z beyond pressure and density measurements



Z is a unique platform for equation-of-state studies



Dynamic material properties (DMP) experiments



- Magnetically launched flyer plates for shock compression¹
 - Flyer impact velocities to ~ 40 km/s
 - Hugoniot states to ~ 10 Mbar; 10,000 50,000 K
 - Pressure and density characterized ~ 1-2 %



¹R.W. Lemke *et al.*, J. Appl. Phys. **98**, 073530 (2005)
²J.-P. Davis *et al.*, Phys. Plasmas **12**, 056310 (2005)
³C. T. Seagle *et al.*, Appl. Phys. Lett. **102**, 244104 (2013)

- Ramp (shockless) compression²
 - Continuous quasi-isentropic compression to ~ 5 Mbar
 - Strain rates ~ 10⁶-10⁷ /s
 - Lower temperature states ~ 1000 3000 K

Shock-ramp compression³

- Initial flyer impact followed ramp loading
- Complex loading path access off-Hugoniot states
- Shock melt and ramp refreeze



Z-DMP planar experiments

- Coaxial load¹
 - Cathode stalk surrounded by anode panels
 - Dual pressures possible on north and south panels
 - Enclosed magnetic fields
 - More sample locations
 - Optimal for (flyer plate) shock compression
- Stripline load²
 - Identical pressure on both cathode and anode panels
 - Higher current density and pressure
 - Open magnetic fields
 - Optimal for high-pressure ramp compression



short circuit

 \vec{B} \otimes

cathode

 $\vec{J} \times \vec{B}$

anode







- ¹M. D. Knudson *et al.*, J. Appl. Phys. **94**, 4420 (2003)
- ²R. W. Lemke *et al.*, Int. J. Impact Eng. **38**, 480 (2011)

Z-DMP cylindrical experiments

- Cylindrical implosion reaches extreme pressure states¹
 - Current pulse shaping creates ramp-wave compression
 - Quasi-isentropic compression to 20 Mbar



$$I = 20 \text{ MA}$$

 $R = 1 \text{ mm}$
 $P_B \approx 64 \text{ Mbar}$

- Diagnostics are challenging²
 - Limited space
 - Miniature probes
 - Velocities well beyond 10 km/s





¹M. R. Martin *et al.*, Phys. Plasmas **19**, 056310 (2012) ²D. H. Dolan *et al.*, Rev. Sci. Instrum. **84**, 055102 (2013)



3 key components to x-ray diffraction on Z-DMP experiments

- Produce source x-rays
 - Laser irradiate metal foil
 - X-pinch
 - X-ray diode
- Generate high-pressure state
 - Z-DMP load
 - Debris mitigation
 - X-ray background
- Detect diffracted x-rays
 - Image plate
 - Scintillator/phosphor
 - Streak camera
 - CCD





Challenges of x-ray diffraction on Z

- Target parameters
 - Large and thick samples
 - Reflection geometry
 - Containment targets
 - Inserting incident x-rays
 - Extracting diffracted x-rays
- Destructive environment of Z-DMP load
 - Prevent catastrophic vacuum breach
 - Protect ZBL
 - Retrieve data
- X-ray background
 - High energy photons (up to 10 MeV) produced
 - Sufficient signal-to-noise
- Electromagnetic pulse (EMP)
 - Fry electronics

Addressing challenges of Z-XRD

- High photon energy (>10 keV), short duration (< 1ns) multi-pulse x-ray source
 - Penetrate into thick targets
 - Temporally resolve phase transformations
- Placing image plate, x-ray CCD, and x-ray streak camera near load
 - Robust x-ray and EMP shielding
 - Advanced debris mitigation
- Convert diffracted x-rays into visible photons
 - X-ray phosphor/scintillator near load
 - Transport light out of load region (fiber or open optics relay)





Laser-driven dynamic x-ray diffraction at Sandia





Exceptional service in the

national

interest

Tommy Ao and Marius Schollmeier

J. Benage, E. Field, M.W. Kimmel, J.L. Porter, P.K. Rambo, J. Schwarz, C.T. Seagle

Sandia National Laboratories, Albuquerque, NM

National ICF Diagnostics Working Group Meeting Los Alamos National Laboratory, Los Alamos, NM, October 6-8, 2015





Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.



Laser and source requirements

Source:

- Above 10 keV due to thick targets
- Monochromatic
- Short emission duration (1 ns or below)
- Multi-pulse with >5 ns inter-pulse delays

 \rightarrow use multiple K α bursts from period-5 transition metals (15-25 keV)

Laser and focusing hardware:

- Multi-pulse capability
- Sub-ns pulse duration
- Final focusing optics well-shielded from Z debris
- Z vacuum protection when debris protection does not hold

 \rightarrow modify ZPW for multi-pulse, 100-ps operation & use existing ZBL lens focusing

Facility Overview





Target Bay Overview



New target chamber (CHAMA)

- 1st chamber to use ZBL, ZPW & CHACO
- anticipated activation: Q1, FY16

Laser systems:

- Chaco laser to load diffraction targets
- Up to 50J/532nm/5-10ns
- Status:
 - Laser operational
 - Beam delivery to CHAMA under construction
- ZPW to create x-ray source
- Up to ~400J/1054nm/50-200ps
- Energy is limited due to gold gratings and Bintegral issues in final focal lens
- Status:
 - Laser operational
 - Optics being coated
 - Mounting hardware under construction





Sub-ns laser-created x-ray sources



- Lens-based focusing of ZPW requires 100-ps-scale pulse duration
 - How efficient is 100-ps x-ray generation?
 - What is the photon yield?
 - How "clean" is the spectrum (He_α)?
 - What is the x-ray pulse duration?
 - What is the optimum target size?

Sandia

National



Debris generation at Z



current ZPW Final Optics Assembly for Z



- ZPW:
 - 43x43 cm beam
 - 500 J, 500 fs
 - F/# = 11 parabola (f = 4.73 m)
- FOA is designed for off-axis, sub-ps irradiation of backlighter targets





J. Schwarz et al., PRST-AB 13, 041001 (2010)

Optics protection for on-axis irradiation





B-Integral considerations

- Refractive index depends on intensity: $n(I) = n_0 + n_2 I$
- Intense laser propagation then has a nonlinear component to phase ($\Phi = 2\pi/\lambda n(l) z$), which accumulates (called B-integral):

$$B = \frac{2\pi}{\lambda} \int_0^L n_2(z) I(z) \, \mathrm{d}z$$

- Limits:
 - Total accumulated nonlinear phase must be ∑B < 4 to avoid whole-beam self-focusing effects (i.e., focal spot shifting)
 - "small" spatial defects can be stripped at pinholes of spatial filters. B-integral between pinholes resets to 0 at each pinhole. Keep ΔB < 2

 \Rightarrow Amplify pulse with 1-1.5 ns chirp, then compress to \approx 100 ps prior to focusing



Expected laser focal spot and intensity

- ZEMAX modeling with the existing ZBL lens to focus a 10 TW (1 kJ/ 100 ps) ZPW beam indicates:
 - $\approx 1.2 \times 10^{10} \text{ W/cm}^2$ at the lens (matches B-Integral model)
 - focal spot: 11 μm × 10 μm (FWHM)
 - ≈ 8 × 10¹⁸ W/cm² at the best focus



Sandia National

horatories

Multi-frame ZPW modification

Sandia National Laboratorie

- Concept:
 - Spatially split beam before injection into main amp section
 - Full delay adjustment (2-20 ns)
 - Separated far-field spots (2 targets)
 - Similar to Z-Beamlet multi-frame backlighter (MFB) concept but without certain energy losses (except for apodization strip at 10% level)
 - Compatible with vacuum compressor in Target Bay
 - Allows 2 beams with each at 500J/0.1ns or 1kJ/0.2ns (B-integral limited)





Multi-frame, time-resolved detector ideas



- First measurements will be time-integrated
 - Use image plate as detector
- later: Scintillator/Phosphor screen + relay optics
 - Single-frame with CCD camera
 - Multi-frame operation with hCMOS camera



- Alternative: DIXI + hCMOS
 - Not clear how close DIXI can be to the target
 - Requires further investigation



END

22



Diagnosing material lattice dynamics during a dynamic compression experiment

- What?
 - Characterize phase transformations that occur in dynamically compressed condensed matter on ns time scales and nm spatial scales
- Why?
 - Such information enables to determinate how material behaves under extreme conditions
 - For most materials, there are very few constraints on existing models for phase transitions under dynamic loading
- How?
 - perform time-resolved, x-ray diffraction measurements on dynamically compressed, polycrystalline matter (dynamic xray diffraction, DXRD)
 - Requires sub-ns x-ray probe to resolve dynamics
 - High-Z material requires >10 keV x-rays
 - Use short-pulse laser to create source
 - Modify short-pulse laser system for dual pulse operation to temporally resolve phase transition dynamics





Z-Backlighter facility laser-to-x-ray conversion efficiency scaling



Sandia National

aboratories

Focusing ZPW with a lens: Multi-pulse operation



Options Considered for ZBL MFB

- Angle multiplexing
 - Via lateral pinhole offset
- Wavelength multiplexing
 Separated by diffraction grating
 - Aperture division multiplexing
 - Separated by wedge refraction
 - Combinations of these



For ZPW MFB, it is best to consider either aperture or angle multiplexing or a combination.



Lens-based focusing



