


Commissioning the P11 Neutron Temporal Diagnostic for High-Neutron-Yield Implosions



J. KATZ, C. STOECKL, J. MAGOON, R. TAYLOR, D. GUY, M. COUCH,
F. EHRNE, D. J. LONOBILE, D. WEINER, E. C. COST, C. REES,
M. ROMANOFSKY, J. SZCZEPANSKI, C. ABBOT, T. LEWIS, and M. MASLYN


University of Rochester, Laboratory for Laser Energetics

Summary


- The P11 neutron temporal diagnostic (NTD) provides precision neutron burn rate and bang-time measurements for high-yield cryogenic implosions
- Port selection at P11 allows for a 9-cm scintillator insertion depth on cryogenic implosions to improve available signal levels and measurement bandwidth
- A 200 \times reduction in neutron noise levels relative to the H5-NTD was achieved by locating the streak camera behind the OMEGA shield wall, 11.4 m from target chamber center (TCC)
- A 16-m optical-path-length image-relay system transports scintillator light from inside the target chamber to the streak camera
- Data is recorded using a Rochester Optical Streak System (ROSS) P5100 streak camera with 15-ps temporal resolution

E23899
UNIVERSITY OF ROCHESTER

The OMEGA Target Bay section includes the scintillator-transport mechanism, zoom-optics assembly, and image-relay hardware




The P11-NTD provides superior data quality on high-yield implosions compared to previous NTD diagnostics




E23906
UNIVERSITY OF ROCHESTER

The NTD measures the neutron burn rate and bang time relative to the OMEGA laser T-0 using a fast-rise-time scintillator




E23900
UNIVERSITY OF ROCHESTER

The final image relay section includes focusing optics, a remote-controlled filter wheel, and the ROSS streak camera




The P11-NTD leverages the power of the ROSS streak-camera platform to provide a well-characterized camera response




E23907
UNIVERSITY OF ROCHESTER

The P11-NTD delivers the instrument performance required to support the current and future LLE cryogenic campaign




E23902
UNIVERSITY OF ROCHESTER

Neutron noise was minimized by placing the streak camera ~11 m from TCC with ~170 cm of line-of-sight shielding



E23905
UNIVERSITY OF ROCHESTER

The P11-NTD maintains excellent CCD image signal-to-noise ratios throughout the typical OMEGA DT neutron yield range




E23919
UNIVERSITY OF ROCHESTER

Summary




- The P11 neutron temporal diagnostic (NTD) provides precision neutron burn rate and bang-time measurements for high-yield cryogenic implosions
- Port selection at P11 allows for a 9-cm scintillator insertion depth on cryogenic implosions to improve available signal levels and measurement bandwidth
- A 200× reduction in neutron noise levels relative to the H5-NTD was achieved by locating the streak camera behind the OMEGA shield wall, 11.4 m from target chamber center (TCC)
- A 16-m optical-path-length image-relay system transports scintillator light from inside the target chamber to the streak camera
- Data is recorded using a Rochester Optical Streak System (ROSS) P5100 streak camera with 15-ps temporal resolution

The NTD measures the neutron burn rate and bang time relative to the OMEGA laser T-0 using a fast-rise-time scintillator




The P11-NTD delivers the instrument performance required to support the current and future LLE cryogenic campaign



Performance metric	Performance status
Minimum burn width	50 ps
Bang-time measurement accuracy	± 50 ps
Detectable DD neutron-yield range	5×10^9 to 1×10^{13}
Detectable DT neutron-yield range	5×10^{10} to 1×10^{15}


E23902

The OMEGA Target Bay section includes the scintillator-transport mechanism, zoom-optics assembly, and image-relay hardware




E23903

The final image relay section includes focusing optics, a remote-controlled filter wheel, and the ROSS streak camera




E23904

Neutron noise was minimized by placing the streak camera ~11 m from TCC with ~170 cm of line-of-sight shielding



E23905

The P11-NTD provides superior data quality on high-yield implosions compared to previous NTD diagnostics



E23906

The P11-NTD leverages the power of the ROSS streak-camera platform to provide a well-characterized camera response



Rochester Optical Streak System
ROSS



Sweep-time base calibration		Camera bandwidth characterization											
Streak-tube geometric distortion	Grid point array 	Cathode image focus time axis 	Static line spread function <p>ADU's ($\times 10^6$)</p> <p>Position (pixel)</p>										
Ramp-driver sweep-speed profile	On-shot 2-GHz comb Dwell time plot <p>Sweep rate (ps/px)</p> <p>CCD position (pixel)</p>	<table border="1"><thead><tr><th>Sweep speed (ns)</th><th>Temporal resolution (ps)</th></tr></thead><tbody><tr><td>13</td><td>60</td></tr><tr><td>6</td><td>25</td></tr><tr><td>3 (nominal)</td><td>15</td></tr><tr><td>1.5</td><td>7</td></tr></tbody></table>		Sweep speed (ns)	Temporal resolution (ps)	13	60	6	25	3 (nominal)	15	1.5	7
Sweep speed (ns)	Temporal resolution (ps)												
13	60												
6	25												
3 (nominal)	15												
1.5	7												
Time base corrected with 1% accuracy		15-ps on-shot camera resolution											

E23907

The P11-NTD maintains excellent CCD image signal-to-noise ratios throughout the typical OMEGA DT neutron yield range



Peak CCD image signal-to-noise (SNR) response

