#### Current Status of Tertiary Neutron Diagnostic by Carbon Activation



Vladimir Glebov et al. University of Rochester Laboratory for Laser Energetics 31st Annual Anomalous Absorption Conference Sedona, AZ 3–8 June 2001



D. D. Meyerhofer, P. B. Radha, W. Seka, S. Skupsky, J. M. Soures, and C. Stoeckl

University of Rochester Laboratory for Laser Energetics

S. Padalino, L. Baumgart, R. Colburn, and J. Fuschino State University of New York at Geneseo

> T. C. Sangster Lawrence Livermore National Laboratory

#### Summary

#### The present carbon activation system can be used for ρR measurements of DT cryogenic targets on OMEGA

- The  $\rho$ R of cryogenic targets can be determined from the tertiary neutron yield.
- We perform several tests on OMEGA to study background and contamination signal of carbon activation.
- For these tests we use DT-filled glass targets with very low  $\rho R$  that don't produce tertiary neutrons.
- Carbon disks were purged in vacuum at 1000°C and filled with different gases: Ar, N<sub>2</sub>, CO<sub>2</sub>, and air.
- Purification and proper packaging reduced the false signal from contamination and secondary reactions to an acceptable level.

#### Outline

# Current status of the tertiary neutron diagnostic by carbon activation



Carbon as an activation material

UR

- Carbon purification
- Experimental setup
- Experimental data
- Interpretation of experimental data
- Conclusion

### Measurement of tertiary neutrons can be used to determine the $\rho$ R of ICF targets

• Primary DT fusion reaction:

 $D + T \rightarrow \alpha + n$  (14.1 MeV)

• Secondary, 14.1-MeV neutrons scatter elastically in the fuel:

 $n + D \rightarrow n' + D' (0-12.5 \text{ MeV})$ 

 $n + T \rightarrow n' + T' (0-10.6 \text{ MeV})$ 

Tertiary, in-flight fusion reaction:

**D'** (0–12.5 MeV) + T  $\rightarrow \alpha$  + n'' (12–30 MeV)

• For NIF, the yield of tertiary neutrons is proportional to  $\rho R$ .

- The (n, 2n) reaction in  $^{12}$ C has a threshold of about 20 MeV.
- <sup>11</sup>C has a half-life of 20.3 min and emits positrons that produce two 511-keV gamma rays upon annihilation.



 Any contamination material that produces a positron emitter by interaction with 14.1-MeV neutrons will be the background for carbon activation.

- Since primary neutron yield is about 10<sup>6</sup> larger than tertiary yield, contamination should be less than one part per million.
- Contamination levels producing the same decay signal as carbon:

Material	Reaction	N <sub>cont.</sub> /N <sub>carbon</sub>
Cu	<sup>63</sup> Cu (n,2n) <sup>62</sup> Cu	1.7 × 10 <sup>-8</sup>
Cu	<sup>65</sup> Cu (n,2n) <sup>64</sup> Cu	3.5 × 10 <sup>−7</sup>
N	<sup>14</sup> N (n,2n) <sup>13</sup> N	8.2 × 10 <sup>-7</sup>
Ni	<sup>58</sup> Ni (n,2n) <sup>57</sup> Ni	1.6 × 10 <sup>−5</sup>
Cr	<sup>50</sup> Cr (n,2n) <sup>49</sup> Cr	1.7 × 10 <sup>−5</sup>

#### A carbon purification facility was built at the State University of New York at Geneseo

- Carbon disks were baked in vacuum at 1000°C for 8 hours.
- The carbon disks were cooled and stayed in a gas atmosphere for 24 hours.
- The carbon disks were then sealed in vacuum bags.
- The vacuum bags were opened 7 to 10 min before the laser shot.



# The OMEGA copper activation system is used for primary DT yield measurements, secondary yield measurements, and carbon activation R & D



Automatic pneumatic retractor Cu/C disk: 76-mm diam  $\times$  9.5 mm Disk-to-target distance: 40 cm

E10676a

3 m

### The present carbon activation system can be used for $\rho R$ measurements of OMEGA cryo-DT targets



Solid angle =  $2.25 \times 10^{-3}$  Activation/incident neutron =  $4 \times 10^{-2}$ Efficiency = 20% Counts/produced neutron =  $1.8 \times 10^{-5}$ At yield =  $10^{14}$  and  $\rho R = 0.1$  g/cm<sup>2</sup> 3600 counts will be detected.

### Carbon activation tests were carried out on the 30-kJ, 60-beam OMEGA laser system

- Glass-shell targets filled with 20 atm of DT were used in implosion experiments.
- These targets have very low  $\rho R$  and do not produce tertiary neutrons.

- A 1-ns square laser pulse shape was used in these studies.
- Primary DT neutron yield was measured by a neutron-time-of-flight (nTOF) scintillating counter.
- The neutron yield in these tests was from  $4 \times 10^{13}$  to  $9.6 \times 10^{13}$ .

# With purification and proper packaging the contamination signal from carbon disk can be reduced substantially



• Expected number of counts from OMEGA DT cryo target is about 2660.

### The importance of "sandwich" can be explained by the process shown here



# Pure samples, packaging, and handling procedures are all important for carbon activation diagnostics

**Different factors in order of importance:** 

- Use pure carbon sample, without internal contamination.
- Carbon disk should be wrapped in a plastic bag to prevent external contamination.
- "Sandwich" carbon disk to shield from protons.
- Minimize disk contact with air before shot.

#### Summary/Conclusion

#### The present carbon activation system can be used for ρR measurements of DT cryogenic targets on OMEGA

- The  $\rho$ R of cryogenic targets can be determined from the tertiary neutron yield.
- We perform several tests on OMEGA to study background and contamination signal of carbon activation.
- For these tests we use DT-filled glass targets with very low  $\rho R$  that don't produce tertiary neutrons.
- Carbon disks were purged in vacuum at 1000°C and filled with different gases: Ar, N<sub>2</sub>, CO<sub>2</sub>, and air.
- Purification and proper packaging reduced the false signal from contamination and secondary reactions to an acceptable level.