#### Direct-Drive Implosion Experiments with Enhanced Beam Balance on the OMEGA Laser



F. J. Marshall, J. A. Delettrez, R. L. Keck, J. H. Kelly, R. Epstein, and L. J. Waxer University of Rochester Laboratory for Laser Energetics 43rd Annual Meeting of the American Physical Society Division of Plasma Physics Long Beach, CA 29 October–2 November 2001 Summary

# A method to enhance the OMEGA laser system's on-target beam balance has been developed

- The on-target beam intensities are determined by analysis of x-ray images of laser-plasma generated by non-overlapping beams.
- The measurements have been used to reduce the beam peak intensity variations from  $\sim 6\%$  (rms) to  $\sim 2\%$  (rms).
- Experiments have been performed with enhanced beam balance, demonstrating improved implosion symmetry.
- Fusion performance (Y<sub>n</sub>, ρR<sub>fuel</sub>) is not significantly improved, however, by enhanced beam balance, implying that other factors such as beam shape or target imperfections may be more important.

#### A method to determine the on-target-beamintensity variation has been developed

- All 60 beams of OMEGA irradiate a 4-mm-diam Au-coated sphere.
- The x-ray emission from each beam spot is corrected for conversion efficiency and view-angle effect to determine the beam-to-beam intensity variations.

 $I_{X} \propto I_{UV}^{\gamma}$  (power-law x-ray conversion dependence)

 $I_{\mathbf{X}}(\theta) = I_{\mathbf{X}}(\mathbf{0})f(\theta)$  (view-angle effect)

 Once the individual-beam-intensity variations are determined, adjustments are made to the last amplifier of each beam to improve on-target balance.

# Each OMEGA pointing shot image is simulated to determine the beam peak intensities



### The x-ray beam spots are fit to super-Gaussian intensity patterns corrected for view angle



### The effect of view angle can be accurately modeled by x-ray emission from an optically thin plasma shell



### Up to eight cross-calibrated x-ray cameras (XPHC's) are used to determine the beam's peak

UR LLE



## The beam-to-beam intensity variations are reduced with enhanced balance

Measured beam peak intensity distributions **Standard balance Enhanced balance** OMEGA shot 24679 OMEGA shot 24769 15 σrms σrms **HED's 2.8% HED's 6.9%** Frequency 10 **XPHC's 6.6% XPHC's 2.2%** 5 0 0.9 0.9 1.2 **8.0** 1.0 1.1 1.2 **0.8** 1.0 1.1 Normalized intensity Normalized intensity

### Enhanced balance implosions obtain more spherically shaped cores

OMEGA direct-drive, D<sub>2</sub>-filled, **18.5-**µm-thick-CH-shell implosions 15-atm-filled 7-atm-filled 3-atm-filled Standard balance  $\leftrightarrow$ **50** µm Enhanced balance **CR** = 14 **CR** = 23 **CR** = 38

# The orientation of perturbations seen in the implosions is consistent with the inferred on-target uniformity

OMEGA direct-drive, D<sub>2</sub>-filled-CH-shell implosions Standard balance (24119) +5% Normalized intensity KB3 **50 μm** Enhanced balance (24121) -5% KB3

#### Target performance is not significantly affected by beam balance with current beam smoothing



Summary/Conclusions

# A method to enhance the OMEGA laser system's on-target beam balance has been developed

- The on-target beam intensities are determined by analysis of x-ray images of laser-plasma generated by non-overlapping beams.
- The measurements have been used to reduce the beam peak intensity variations from ~6% (rms) to ~2% (rms).
- Experiments have been performed with enhanced beam balance, demonstrating improved implosion symmetry.
- Fusion performance (Y<sub>n</sub>, ρR<sub>fuel</sub>) is not significantly improved, however, by enhanced beam balance, implying that other factors such as beam shape or target imperfections may be more important.