


Neutron Burn History Measurements of D₂ Cryogenic Targets on OMEGA



V. Yu. Glebov
University of Rochester
Laboratory for Laser Energetics

44th Annual Meeting of the
American Physical Society
Division of Plasma Physics
Orlando, FL
11–15 November 2002

Collaborators



**J. A. Delettrez, P. W. McKenty,
P. B. Radha, T. C. Sangster, and C. Stoeckl**

**Laboratory for Laser Energetics
University of Rochester**

Related Talks:

**J. A. Delettrez GO2.004
T. C. Sangster RI1.006**


Summary

CryoNTD records the important neutron burn history information for D₂ cryogenic implosions on OMEGA




- CryoNTD was built because the NTD system on OMEGA is too insensitive to measure the burn history of D₂ cryogenic implosions due to standoff requirements of the Cryogenic Target Handling System.
- CryoNTD has a time resolution ~ 80 ps, absolute timing calibration ~ 40 ps, and a sensitivity of $\sim 10^9$ neutrons.
- CryoNTD data is correlated with ice-layer quality (melting or crystallization) at shot time.
- The neutron burn history of the layered D₂ cryogenic targets generally agrees with the *LILAC* calculation.

We built a TIM-based cryoNTD specially for D₂ cryogenic target shots on OMEGA




NTD
0.6-cm scintillator
at 20 cm from TCC




CryoNTD
3-cm scintillator
at 9 cm from TCC

The cryoNTD setup uses the OMEGA fiducial system as a timing reference




Cryo NTD sensitivity is about 1×10^9 neutrons


Streak camera image



Shot 27143
 $Y = 1.72 \times 10^9$




Absolute timing of the cryoNTD was established with 40-ps accuracy using NTD as a reference



A layered cryogenic target can be in one of the three stages inside the OMEGA target chamber


A cryogenic target inside the target chamber is in the layering sphere before the shot.



Heating = Cooling
layered target



Heating < Cooling
"Crystallized" target




Heating > Cooling
"Melted" target




Transmission of the IR fiber can change up to 10% during insertion into the target chamber.


A “crystallized” cryogenic target has a neutron burn width much wider than predicted by *LILAC*




A “melted” cryogenic target has a neutron bang time much earlier than predicted by *LILAC*



A layered cryogenic target has a neutron burn history in general agreement with *LILAC* prediction



Layered cryogenic target has neutron burn history in general agreement with *LILAC* prediction



Summary/Conclusions

CryoNTD records the important neutron burn history information for D₂ cryogenic implosions on OMEGA



- CryoNTD was built because the NTD system on OMEGA is too insensitive to measure the burn history of D₂ cryogenic implosions due to standoff requirements of the Cryogenic Target Handling System.
- CryoNTD has a time resolution ~ 80 ps, absolute timing calibration ~ 40 ps, and a sensitivity of ~ 10⁹ neutrons.
- CryoNTD data is correlated with ice-layer quality (melting or crystallization) at shot time.
- The neutron burn history of the layered D₂ cryogenic targets generally agrees with the *LILAC* calculation.