


Cryogenic Target Characterization at LLE



D. D. Meyerhofer
University of Rochester
Laboratory for Laser Energetics

45th Annual Meeting of the
American Physical Society
Division of Plasma Physics
Albuquerque, NM
27–31 October 2003

Collaborators



**W. Seka, M. Alexander, R. S. Craxton, M. D. Wittman,
M. Pandina, L. S. Iwan, L. M. Elasky, D. R. Harding, S. Jin,
T. J. Kessler, R. L. Keck, L. D. Lund, and D. Weiner**

**University of Rochester
Laboratory for Laser Energetics**

A. Warrick

**Lawrence Livermore National Laboratory
(was on assignment at LLE)**

T. Brown

**University of Rochester
The Institute of Optics**

Cris Cotton

Optics Consultant


Summary

LLE has made significant progress in characterizing cryogenic targets



- 1-D ray tracing and 3-D optics codes are used to understand the imaging.
- Interpretation routines can handle spider web suspension and other surface perturbations.
- Multiple views during layering provide
 - average power spectra
 - 3-D surface images
 - spherical-harmonic decomposition ($Y_{\ell m}$)
- Targets are also characterized immediately prior to the shot.

LLE cryogenic target characterization is based on a diffuse f/15 source and f/6 (f/15) imaging optics




Cryogenic targets are located inside the cryogenic shroud with three windows at input and three at output ($\sim f/6$ limit).

A pulsed LED light source and a high-quality CCD camera are used for target characterization



- High-brightness LED at 644 nm \square 0.2 ms exposures
- DALSA CCD triggers on LED pulse for timing and integration
- Darkfield and flatfield corrections improve image quality
- Rapid image transfer CCD permits images to be taken during last 50 ms before shot with minimum interference at shot time
- Sapphire targets serve to determine lowest order optical image distortions \square appropriate corrections are then applied to all images
- Sapphire targets also serve for calibration of target sized and measured target distortions via AFM-characterization



Various rings focus at different locations and encode different information that has been identified with the 1-D ray trace code



- In real cryogenic targets, the bright ring is located further inside, and the other rings are more closely spaced but still visible.


Recent D₂-ice layers with IR heating approach direct-drive ICF requirements at 1 to 2 times the equivalent tritium decay heating (Q_{DT})

Average spectra for 23 different angular views of the same cryo target




930- μm -diam OMEGA cryo target with 100- μm -D₂-ice layer and 3.5- μm -CH shell

Different views give different ice roughnesses and spectra



The roughnesses from shadowgraphy are larger than from AFM – under investigation


Three-dimensional reconstructions have been obtained for cryogenic targets




- Surface position is mapped onto sphere.
- Data are smoothed.
- Information for low-order modes is provided.

A spherical-harmonic ($Y_{\ell m}$) decomposition up to $\ell = 8$ and including all m 's has been carried over

3-D reconstruction




Result of $Y_{\ell m}$ decomposition




A spherical-harmonic ($Y_{\ell m}$) decomposition up to $\ell = 8$ and including all m 's has been carried over

Result of $Y_{\ell m}$ decomposition




$A_{\ell m}$ coefficients for $\ell_{\max} = 8$,
target: cryo2038-220



Shadowgraphs are also taken of targets at shot time



- A shadowgraphy system has been installed on the OMEGA target chamber.
- The target is observed after shroud removal, ~10 ms before the shot.
- Even though it is a static camera, the finite readout allows vibrational amplitude to be inferred.



LLE has made significant progress in characterizing cryogenic targets

- 1-D ray tracing and 3-D optics codes are used to understand the imaging.
- Interpretation routines can handle spider web suspension and other surface perturbations.
- Multiple views during layering provide
 - average power spectra
 - 3-D surface images
 - spherical-harmonic decomposition ($Y_{\ell m}$)
- Targets are also characterized immediately prior to the shot.