


Update on the Rochester Optical Streak System



**P. A. Jaanimagi, R. Boni,
and D. D. Meyerhofer
University of Rochester
Laboratory for Laser Energetics**

**47th Annual Meeting of the
American Physical Society
Division of Plasma Physics
Denver, CO
24–28 October 2005**

Summary

The ROSS camera is a comprehensive diagnostic system with autofocusing and self-calibration capability




- The high-precision measurements demanded by the LLE experimental program require precise setup and calibration of the diagnostics.
- The optical calibration module incorporates an extensive imaging, flat-field correction, geometric-distortion, and time-calibration capability.
- All functions can be accessed and monitored remotely.

All pixels in the recorded image were not created equal




- **The system response is dependent on the spatial and temporal position and the focusing conditions.**
- **To make pretty pictures and FWHM measurements requires minimal calibration.**
- **Serious quantitative measurements (SNR > 10, DR > 10, nonlinearity < 10%) require extensive calibrations.**
- **These issues are addressed in the ROSS design with the optical calibration module (OCM).**

The input imaging system is an Offner triplet with motorized controls for the dual object planes and the secondary mirror



Flat-fielding a CCD recording system elevates its performance to a near-quantum-limited detector




The ROSS–OCM design uses an LED, homogenizer, and a slow ramp to flat-field the system



- Input optics, photocathode quantum efficiency, photoelectron transport, phosphor-conversion efficiency, fiber-optic coupling, and CCD quantum efficiency are all independent of the signal level.
- Uses a 1- to 10-s duration ramp to simplify the illumination-source design.
- Multiple-wavelength CW–LED's are available internally, or white light passed through an interference filter can be fed in via fiber optics.
- Homogenizer output is spatially uniform to 0.25%.


Geometric distortions result from electron-optic aberrations and mechanical misalignments

- Manufacturers of streak tubes are reluctant to address these issues—\$\$
- Use ROSS–OCM flat-field subsystem plus spatial and temporal modulation to produce raw uncorrected images
- $S(x,y) = \text{Jacobian}(x,y,x',y') \times S(x',y')$, $x = U(x',y')$ and $y = (x',y')$



Time calibration is accomplished with an on-board, 664-nm, 3-GHz comb generator

- Verify system operation prior to shot
- External source may also be coupled into any or all channels




**4 streak speeds
(2–25 ns) can be
selected remotely**


The autofocus capability allows local optimization of the spatial or temporal

- Geometric distortion image shows the main EO aberration is field curvature.

Annulus is the intersection of a curved image plane with a flat screen




Lineout through central region



The P820 streak tube provides a 2-ps time resolution capability for OMEGA EP diagnostics

- The peak current handling capability is 1.8 mA.



Summary/Conclusions

The ROSS camera is a comprehensive diagnostic system with autofocusing and self-calibration capability



- The high-precision measurements demanded by the LLE experimental program require precise setup and calibration of the diagnostics.
- The optical calibration module incorporates an extensive imaging, flat-field correction, geometric-distortion, and time-calibration capability.
- All functions can be accessed and monitored remotely.