#### First Tests of a Bubble Chamber for Neutron Imaging on OMEGA



M. C. Ghilea University of Rochester Laboratory for Laser Energetics 49th Annual Meeting of the American Physical Society Division of Plasma Physics Orlando, FL 12–16 November 2007

Summary

### A liquid-freon neutron detector is being tested for use in a neutron imaging system on OMEGA

- The detection system has been assembled.
- Individual bubbles have been imaged to study bubble growth.
- The published\* phase behavior of freon has been confirmed.
- The operational window for producing bubble images has been determined.

#### **Collaborators**



D. D. Meyerhofer, T. C. Sangster, D. J. Lonobile, and A. Dillenbeck Laboratory for Laser Energetics University of Rochester

R.A. Lerche

Lawrence Livermore National Laboratory

L. Disdier

Commissariat à L'Énergie Atomique

### A bubble chamber will be placed in the target bay 7.6 m from TCC to measure the neutron source distribution

### Initial tests will determine • Resolution, using hard edges at 7.6 m the detector • Flat-field response TIM 6 26.5°

# The encoded neutron distribution from an aperture will be measured by the bubble density distribution

- After decompressing liquid freon below the liquid equilibrium line, a neutron interaction creates a bubble.
- The bubble location gives the same information as a pixelated scintillator array (i.e. spatial location).
- The bubble density is expected to be too high for individual imaging, so a Schlieren method will be used.

## The LLE bubble chamber has been designed to measure the column density of bubbles



#### The operation pressure cycle depends on the phase behavior of freon



The detector operates in cycles between the "start" and "working" points.

### The bubble growth rate is determined by the freon temperature



UR

The temperature is tuned to optimize the bubble size given readout duration.

150

100

50

0

Bubble diameter (µm)

# Test images show the importance of the time sequence of the decompression cycle and the probe laser with respect to $T_0$



The rising bubbles leave a thermal trail that will likely affect the image resolution.

## There is a window of about 40 ms for imaging before the freon begins to boil



The probe laser is pulsed for 1 ms and the relative accuracy of timing with OMEGA can be controlled to within 0.1 ms.

Summary/Conclusions

#### A liquid-freon neutron detector is being tested for use in a neutron imaging system on OMEGA

- The detection system has been assembled.
- Individual bubbles have been imaged to study bubble growth.
- The published\* phase behavior of freon has been confirmed.
- The operational window for producing bubble images has been determined.