



High-Resolving-Power, Ultrafast Streaked X-Ray Spectroscopy on OMEGA EP



P. M. Nilson
University of Rochester
Laboratory for Laser Energetics

58th Annual Meeting of the
American Physical Society
Division of Plasma Physics
San Jose, CA
31 October–4 November 2016

A high-resolving-power, streaked x-ray spectrometer is being developed and tested on OMEGA EP



- The goal is to achieve a resolving power of several thousand and 2-ps temporal resolution (February 2017)
- To understand system performance, a time-integrating survey spectrometer has been deployed on OMEGA EP
- Survey spectrometer measurements and offline testing show
 - focusing fidelity: $\sim 50\text{-}\mu\text{m}$ line focus
 - resolving power: >2000
 - throughput: $\sim 10^{-7} \text{ ph/ph}$
 - shielding: 5 to 15 cm of lead
- These measurements provide a firm foundation for designing and implementing the time-resolved instrument

The instrument will ultimately be used to measure temperature equilibration dynamics and material response to ultrafast heating.

Collaborators



**F. Ehrne, C. Mileham, D. Mastrosimone, R. K. Jungquist,
C. Taylor, R. Boni, J. Hassett, C. R. Stillman, S. T. Ivancic, D. J. Lonobile,
R. W. Kidder, M. J. Shoup III, A. A. Solodov, C. Stoeckl, and D. H. Froula***

**University of Rochester
Laboratory for Laser Energetics
*also Department of Physics**


**K. W. Hill, L. Gao, M. Bitter, and P. Efthimion
Princeton Plasma Physics Laboratory**

**D. D. Meyerhofer
Los Alamos National Laboratory**

A high-energy ultrafast laser can heat solid-density material on a time scale much faster than the material expands

- Heating at high density produces exotic states of matter in extreme thermodynamic conditions¹
- The possible extremes in temperature enables novel material and radiative properties experiments²
 - e.g., mean opacity of solar interior matter³
- New diagnostic techniques are sought for testing
 - plasma-dependent atomic processes⁴
 - plasma opacity⁵
 - equation-of-state models⁶

These studies require dense, high-temperature plasmas that are well characterized.



¹ A Report on the SAUUL Workshop, Washington, DC (17–19 June 2002).

² K. Nazir et al., Appl. Phys. Lett. **69**, 3686 (1996).

³ J. E. Bailey et al., Nature **517**, 56 (2015).



⁴ D. J. Hoarty et al., Phys. Rev. Lett. **110**, 265003 (2013).

⁵ R. A. London and J. I. Castor, High Energy Density Phys. **9**, 725 (2013).

⁶ M. E. Foord, D. B. Reisman, and P. T. Springer, Rev. Sci. Instrum. **75**, 2586 (2004).

Outer-shell ionization affects the energy and shape of the characteristic K_{α} line in a partially ionized plasma


- Hot electrons create K-shell vacancies when colliding with ions
- Ionization by thermal electrons removes electrons from the ions' outer shells
- As the ionization progresses, the $K_{\alpha_{1,2}}$ lines increase their energy^{1–5}



Time-resolving the K_{α} line shift allows for the mean ionization state of the plasma to be inferred during the rapid heating phase.


- ¹K. Słabkowska et al., High Energy Density Phys. **15**, 8 (2015).
²K. Słabkowska et al., High Energy Density Phys. **14**, 30 (2015).
³G. Gregori et al., Contrib. Plasma Physics **45**, 284 (2005).
⁴P. M. Nilson et al., Phys. Plasmas **18**, 042702 (2011).
⁵J. F. Seely et al., High Energy Density Phys. **9**, 354 (2013).

The instrument is based on two diagnostic channels, each with a spherical Bragg crystal




Parameter	Requirements
X-ray source size	$\sim 100 \mu\text{m}^2$
Spectral range	7.97 to 8.11 keV
Crystal and Bragg angle	Si220 crystal—Bragg angle = 22.8°
Crystal radius of curvature	330 mm
Crystal size	25 mm × 100 mm
Source-to-crystal distance	2.2 m
Resolving power	~5000—streak-camera limited
Spectral shifts	Few eV to 20-eV K_{α} line shifts
Streak-camera slit	6-mm-long, 400- μm -wide 50- μm -high-throughput region
Temporal resolution	2 ps

High-power experiments show the focusing fidelity, resolving power, and throughput meet the desired requirements



Phase I has provided the foundation for designing and implementing the time-resolved instrument.

Time-integrated measurements on OMEGA EP show spectral shifts increasing with target energy density



A high-resolving-power, streaked x-ray spectrometer is being developed and tested on OMEGA EP




- The goal is to achieve a resolving power of several thousand and 2-ps temporal resolution (February 2017)
- To understand system performance, a time-integrating survey spectrometer has been deployed on OMEGA EP
- Survey spectrometer measurements and offline testing show
 - focusing fidelity: $\sim 50\text{-}\mu\text{m}$ line focus
 - resolving power: >2000
 - throughput: $\sim 10^{-7} \text{ ph/ph}$
 - shielding: 5 to 15 cm of lead
- These measurements provide a firm foundation for designing and implementing the time-resolved instrument

The instrument will ultimately be used to measure temperature equilibration dynamics and material response to ultrafast heating.

Temporal spectral shifts on the Cu K α line in rapidly heated solid matter will validate the spectrometer performance

- Synthetic spectra from hot, dense matter are required
- *LSP*¹ calculates
 - energy-transport physics
 - electromagnetic-field generation
 - target heating
- *LSP* is post-processed based on tabulated *PrismSPECT*² calculations using
 - the local density and temperature at the time of emission
 - line-of-sight and high- T_e opacity effects
- The calculations use an occupation probability model³ and the ionization potential depression formalism of More⁴



¹D. R. Welch et al., Phys. Plasmas **13**, 063105 (2006).

²Prism Computational Sciences Inc., Madison, WI 53711.

³D. G. Hummer and D. Mihalas, Astrophys. J. **331**, 794 (1988).

⁴R. M. More, J. Quant. Spectrosc. Radiat. Transf. **27**, 345 (1982).