Broadband Reflectivity Diagnostic Development for Dynamic Compression Experiments on OMEGA EP



B. Henderson University of Rochester Laboratory for Laser Energetics 60th Annual Meeting of the American Physical Society Division of Plasma Physics Portland, OR 5 - 9 November, 2018

We are developing a white light reflectivity diagnostic to explore the electronic properties of matter at extreme conditions

- HED reflectivity measurements are mostly monochromatic (VISAR), requiring assumptions about the wavelength dependence of electronic behavior
- We have developed a white light source by using a 9-eV shocked silica aerogel
- The reflected white light probe is used to measure the broadband (400 to 800nm) reflectivity of ramp-compressed materials
- Preliminary experiments show sufficient signal-to-noise to resolve changes in reflectivity of matter at extreme conditions





T. R. Boehly, M. Zaghoo, J. R. Rygg, D. N. Polsin, X. Gong, L. Crandall, M. Huff, M. K. Ginnane, and G. W. Collins

> University of Rochester Laboratory for Laser Energetics

> > S. Ali and P. M. Celliers

Lawrence Livermore National Laboratory



Motivation

We will use a broadband reflectivity diagnostic to investigate the electronic properties of HED matter



Backlighter Design

A silica aerogel backlighter is shocked to ~9 eV, producing a bright white light source



Shocked silica foam uses a less energetic input to maintain a steady, high-temperature (9-eV) profile over several nanoseconds



An *in-situ* optical train directs the white light probe onto the back surface of the sample

- The backlighter's shock emission is directed to the back surface of the sample
- The reflected signal is collected with an achromatic *f*/3.3 telescope and relayed to a steak camera
- VISAR simultaneously records the sample's particle velocity





Method

A reflective reference layer is embedded in the sample stack to provide a constant reflectivity baseline for the measurement





Method

The reflected signal is filtered with a stack of narrowband filters and streaked in time





Results

The SOP streak image shows the time-dependent reflected signal from the sample and reference



Results

A simultaneous measurement with VISAR allows us to determine the pressure in the sample



• An iterative Lagrangian analysis (ILA) determines the pressure in the Sample-LiF interface



LLE

Results

The signal-to-noise of this method is sufficient to resolve time-dependent changes in reflectivity





Summary/Conclusion

We are developing a white light reflectivity diagnostic to explore the electronic properties of matter at extreme conditions

 HED reflectivity measurements are mostly monochromatic (VISAR), requiring assumptions about the wavelength dependence of electronic behavior

- We have developed a white light source by using a 9-eV shocked silica aerogel
- The reflected white light probe is used to measure the broadband (400 to 800nm) reflectivity of ramp-compressed materials
- Preliminary experiments show sufficient signal-to-noise to resolve changes in reflectivity of matter at extreme conditions
- Future experiments will incorporate a spectrometer for continuous reflectivity measurements in the wavelength domain



